

ASXL1 p.G646Wfs*12 Detection Enabled and Integrated into a Rapid Automated Myeloid Variant Profiling Solution on Ion Torrent™ NGS Genexus™ System

Brendan Deal, Cheng-Zong Bai, Jingwei Ni, Jennifer Burke, Sarah Brozio, Nitisha Thapa, Jiajie Huang, Sophie Rozenzhak, Seth Sadis, and Marina Sedova, Thermo Fisher Scientific, 200 Oyster Point Blvd, South San Francisco, CA, 94080

INTRODUCTION

Comprehensive gene variant profiling using next-generation sequencing (NGS) improved the research of relevant biomarkers for myeloid disorders. Previously, we described the capability of the Genexus™ System to deliver high quality Variant Calling (VC) results with the OncoPrint™ Myeloid Assay GX v2 (OMAv2) in a 24h automated sample-to-report workflow with minimal hands-on-time.⁽¹⁾

However, certain variants in low complexity genomic regions, like homopolymers, may pose challenges for sequencing and analysis, requiring additional testing and delaying delivery of variant results. One such mutation commonly observed in myeloid malignancies is ASXL1 p.G646Wfs*12 c.1934dupG COSM34210 (Figure 1). Therefore, we explored general improvements to the GX OMAv2 including the ability to detect ASXL1 c.1934dupG.

Figure 1. Schematic of ASXL1 p.G646Wfs*12 c.1934dupG variant

A. ASXL1 wild-type (8G c.1927-1934)

5' ATC GGA GGG GGG GGT GGC CCG 3'
I G G G G G P

B. ASXL1 c.1934dupG (9G)

5' ATC GGA GGG GGG GGG TGG CCC 3'
I G G G G W P

The ASXL1 c.1934dupG variant results from the insertion of an additional guanine next to an 8G homopolymer in the wild-type sequence (A). This insertion causes a frameshift and converts a glycine amino acid to a tryptophan (B). The variant is often challenging to detect due to the low complexity of this region of the gene.

MATERIALS AND METHODS

250+ BAM files were obtained with sequencing DNA from healthy blood donors, commercially available DNA analytical controls, and DNA research samples with known presence of ASXL1 c.1934dupG. The dataset was used for fine-tuning of variant specific parameters in Torrent Variant Caller (TVC) algorithm to achieve >95% Sensitivity and >95% Specificity for ASXL1 c.1934dupG at variant allele frequency (VAF) >10%.

To evaluate the effect of the new TVC parameters, we tested two Myeloid DNA analytical controls, Horizon Myeloid gDNA reference standard HD829 and Seraseq™ Myeloid DNA Mutation Mix (SCMM), containing ASXL1 c.1934dupG, in multiple replicates (N≥20). The testing was performed in nucleic acid-to-report configuration and Genexus™ Software output was used for VC analysis.

HD829 was mixed with variant negative samples, NA12878 or NA24385, to target ASXL1 c.1934dupG VAF below 20%. DNA mixtures were sequenced with ≥5 replicates per VAF in the range 5-18%.

ASXL1 c.1934dupG variant detection was also evaluated in Acute Myeloid Leukemia (AML) research blood samples, herein referred as AML-1, AML-2, and AML-3. Sanger sequencing was used as an orthogonal method to confirm the variant presence. OMAv2 gene specific primers targeting amplicon with ASXL1 c.1934dupG variant were 5' and 3' tailed with an M13 universal sequences to make them compatible with Applied Biosystems™ Dye Primer sequencing. VAFs from the Sanger sequencing data were estimated using QSVanalyser.

Key ASXL1 variant in low complexity region detected by OncoPrint™ Myeloid Assay GX v2

RESULTS

Figure 2. ASXL1 c.1934dupG detection in Myeloid Controls (A) and in AML samples (B)




Figure 3. Correlation of reported and targeted VAF




Figure 5. Sanger sequencing confirms ASXL1 c.1934dupG presence in AML samples but lacks sensitivity to call variant in Seraseq Myeloid DNA Mutation Mix




Table 1. VAF in AML samples estimated with QSVanalyser and reported by Genexus

Sample	Forward strand (%)	Reverse strand (%)	Reported by Genexus
AML-1	37.8	23.8	20.3 (n=6)
AML-2	39.1	20.2	19.1 (n=11)
AML-3	40.4	25.1	25.3 (n=11)
SCMM	7.1	NA	6.5 (n=24)

- VAF values in Sanger sequencing results were estimated to be between 37-40% in the forward DNA strand and 20-25% on the reverse strand.
- For SCMM VAF was estimated to be ~7% with DNA forward strand. DNA reverse strand sequencing data were too noisy to estimate VAF.
- Note the systematic difference between forward and reverse strand VAF.

False Positive (FP) rate and Specificity evaluation

FP rate was estimated using sequencing results from 550+ variant negative samples. Five FP were reported yielding FP rate below 1% and Specificity >99%. The average FP VAF was 4.4%. Yannakou, C.K., Jones, K., McBean, M. et al. ⁽²⁾ discuss the limitation of different detection technologies, including Sanger sequencing, fragment analysis, and NGS.

CONCLUSIONS

Variant calling improvements were made in the OncoPrint™ Myeloid Assay GX v2 to enable robust detection of challenging variants in low complexity regions including ASXL1 c.1934dupG. The capability to accurately report ASXL1 c.1934dupG within OMAv2 is important to provide complete results. Our data support an LOD of 15% VAF for the detection of this variant with the improved analysis workflow. Using the Genexus™ automated NGS sequencing system, the OMAv2 assay enables detection of key Myeloid variants with a 24hr turn-around-time.

REFERENCES

- Marina Sedova, Alexy Ongpin, Jennifer Burke, Collyn Seeger, Sarah Brozio, Janice Au-Young, Jiajie Huang, Thilanka Jayaweera, Iris Casuga, Milton Huynh, Fiona Hyland. Fully automated sample-to-report NGS workflow for comprehensive genomic profiling for myeloid neoplasms [abstract]. In: Proceedings of the American Association for Cancer Research Annual Meeting 2022; 2022 Apr 8-13. Philadelphia (PA): AACR; Cancer Res 2022;82(12_Suppl):Abstract nr 764.
- Yannakou, C.K., Jones, K., McBean, M. et al. ASXL1 c.1934dup;p.Gly646Trpfs*12—a true somatic alteration requiring a new approach. Blood Cancer Journal 7, 656 (2017)

ACKNOWLEDGEMENTS

Thank you to all the individuals involved in the development, commercialization, and constant improvement of the OncoPrint™ Myeloid Assay GX v2.

COPYRIGHT/TRADENAME DISCALIMER

© 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. SeraSeq is a trademark of SeraCare Life Sciences, Inc.